

Freeze Crusher μ T-48

Powerful crushing of frozen samples with liquid nitrogen. 48 samples can be treated simultaneously. Optimum for extracting hard samples, proteins susceptible to heat denaturation, RNA, etc.

Example of various frozen crushed samples including inanimate objects. --> P.091-092



Features

- Crushing of frozen sample in vessels with liquid nitrogen
- 2mL Microtube or Dedicated metal container are used
- The throughput is 0.2g to 2g (Depends on the vessels)

Applications

- Crushing of Yeast, Mold, Tissue piece of animals and plants
- Crushing of bones, teeth and limbs of small animals
- Crushing of Wire covering and Plastics, Asbestos sample etc.

The procedure of freeze crushing Work gloves are worn when freezing the sample with liquid nitrogen in the photo but use gloves suitable for handling liquid nitrogen in actual use.



Put the specified amount (below) sample and metal crusher into the vessel.




When a vessel requires holder, attach it in the holder.



Soak and freeze in liquid nitrogen until bubbling ceases.



Quickly attach to the μ T-48 main unit.



Close the hood, Check the time and shaking speed and Start crushing.



After crushing, take out the crusher and proceed to the next process.

Adapted to various samples with 3 types vessel holders

Adapted to various samples with 3 types vessel holders. Powerful crushing of frozen samples with liquid nitrogen. Living tissues and organs, Hard tissues such as bones etc. and also some inanimate samples such as rubber and plastic can be crushed. Adapted to Marketed 2mL tubes and Dedicated stainless steel crushing vessels.

Optional accessories: Vessel holders



Product Name / Model / Remarks

①48pcs-Holder for μ T-48 TH-0248T

1pc of Holder (Capacity : 48pcs of 2.0mL round bottom Microtubes) and 100pcs of Metal crusher come as a set.

②3pcs-Holder for μ T-48 TH-0203T

4pcs of Holder (Capacity : 3pcs of 2.0mL round bottom Microtubes), 24pcs of Metal crusher and Rack come as a set.

③Stainless steel-made strong crushing vessel TH-SPT

Crushing vessel 4pcs, Dedicated crusher and Rack come as a set. Suitable for samples that cannot be crushed by Microtube with Metal crusher. Larger amount of samples can be crushed than that of Microtube.

Recommended shaking speed in each crushing sample.

- Stainless steel strong crushing vessel: Up to 1000r/min
- Metal crusher: Up to 1200 r/min
- Crushing beads: Up to 1600 r/min

Shaking speed more other than above speed might cause breakage of tubes, vessels. Thus, be sure to observe the shaking speed above.

USER'S VOICE

Very useful when extracting substances susceptible to denaturation and degradation of RNA and proteins.



Model	μ T-48
Crushing method	Vertically reciprocal shaking.
Shaking speed (*1)	0 to 2500 r/min (*1)
Capacity	2.0 mL Microtube: Max. 48pcs (*2)(*3) Stainless steel powerful crushing vessel : Max. 4pcs (*3)
Timer	1 to 999 seconds
Safety devices/ functions	Holder attachment detection switch, Cover opening detection switch
Dimensions	220(W) x 310(D) x 315(H)mm
Weight	Approx. 10.0kg
Power supply	AC100V/1A

(*1) Around 1200r/min should be necessary and sufficient condition to crush the sample in actual use.

(*2) Eppendorf "Safe-Lock Tube 2mL" is recommended.

(*3) Microtube and Stainless powerful crushing vessel are available as an option.

Example ① Freeze crushing of various samples including inanimate samples.

Embrittlement by freezing enables strong crushing.
The freeze crushing with μ T-48 also suitable for Obligatory anaerobe samples.




Freeze Crusher μ T-48 --> P.090

Test results

- Freezing method Immerse the vessels with sample and crusher into liquid nitrogen (2.0 mL: Vessel holder) and freeze them.
- Rushing time 30 sec (Further 30 sec if not completely crushed).
- Judgment whether sample crushed. Whether powder form or nearly it (Cut samples to any size that can be put in vessels).
- Vessels Safe-Lock tube 2.0mL.....Marketed product(Eppendorf made)
 Metal crusher.....Included in Optional 48pcs-holder for μ T-48 (used in this experiment).
 Stainless steel-made strong crush vessel VesselsOptional parts(Dedicated crusher comes with).

Sample	Crushed State	Vessels	Sample volume	Shaking speed	Crushed with
Chicken thigh		Safe-Lock tube 2.0 mL	0.1 g	1200 r/min	Metal crusher
Human hair		Safe-Lock tube 2.0 mL	0.1 g	1200 r/min	Metal crusher
Human nails		Safe-Lock tube 2.0 mL	0.2 g	1200 r/min	Metal crusher
Mouse skin (with body hair)		Safe-Lock tube 2.0 mL	0.2 g	1200 r/min	Metal crusher
Mouse heart		Safe-Lock tube 2.0 mL	0.2 g	1200 r/min	Metal crusher
Mouse tail		Stainless steel-made strong crush vessel	1g	1000r/min	Dedicated crusher
Hypocotyl of Radish		Safe-Lock tube 2.0 mL	0.2 g	1200 r/min	Metal crusher
Okra seeds		Safe-Lock tube 2.0 mL	2	1200 r/min	Metal crusher
Hard rubber (Polychloroprene)		Stainless steel-made strong crush vessel	2g	1000r/min	Dedicated crusher

Frozen sample/Holder temperature (Reference)



● Protuberances not included in Dimensions. ● Vessels of photo not included.

- Selection guide
- Constant temperature incubator shaker OD-Monitor
- CO₂ incubator shaker
- Shaker
- Mixer Rotator Stirrer
- Bead beater homogenizer Ultrasonic homogenizer
- Aluminum block bath Minimize bath
- Water bath Shaking water bath Immersion cooler
- Hybridization oven Constant temperature chamber
- Centrifugal concentrator Cold trap
- Freeze dryer
- Electrophoresis and Blotting apparatus
- Constant temperature water circulating system [Chiller]
- Appendix

Example ② Freeze crushing of various samples including inanimate samples

Freeze crushing of Plastic samples using Freeze crusher μ T-48 with Stainless steel-made strong crush vessel.

Freeze Crusher μ T-48 --> P.090



Results and Examination

We tried some crushing of samples such as polystyrene, polypropylene and polycarbonate that were well known. Each result is as follows.

Sample	Result	Sample shape and volume	Shaking speed	Crushing time
Polystyrene ① : Fair		φ 6 mm ball shaped, 1g	1000 r/min	300 sec
Polystyrene ② : Fair		φ 6 mm ball shaped, 0.4g	1000 r/min	300 sec
Polystyrene ③ : Excellent		10 mm square chip shaped, 0.5 g	1000 r/min	180 sec
Polypropylene ① : Good		10 mm square chip shaped, 0.5g	1000 r/min	150 sec
Polypropylene ② : Good ?		10 mm square chip shaped, 0.5 g	1100 /min	300 sec
Polycarbonate : Poor		10 mm square chip shaped, 0.5 g	1000r /min	300 sec

Polystyrene were able to be powdered completely (③). However, in case of φ6mm ball- shaped sample, there were large fragments remained at certain rate even after trying with different amount and crushing time (①②). It seems that ball-shaped sample remain uncrushed if it is stuck in upward of the crusher.

The result suggests that the shape of sample is better to be chip like shape(or tablet like shape) Polypropylene was crushed into fine fragments although it was not crushed into powder (①).

In order to improve (although the shaking speed limit is exceeded when using a strong crushing container), when it is performed at 1100 r / min for 5 minutes, it becomes fine but braided piece of cotton (②). At this stage, it becomes difficult to collect unless suspended in a solvent. Polycarbonate proved to be difficult to crush. Even if the shaking speed was reduced to 1100 r / min or reducing sample amount, the result of this experiment was that only a small amount of powder was produced and the chip shape remained almost unchanged.

Freeze crushing procedure when using stainless steel-made strong crushing vessel

Adjustable plier is useful for taking in and out of Stainless steel-made strong crushing vessel (referred to as crushing vessel) into liquid nitrogen. Be sure wear leather gloves when touching frozen crushing vessel or the vessel holder that has become cold by contacting the frozen vessel. Be sure ventilate the room well when using liquid nitrogen as there is risk to get Anoxia unknowingly because the vaporized liquid nitrogen becomes huge volume of nitrogen gas.

- Pour liquid nitrogen into a styrofoam container. *1
- Put the weighed sample * 2 into the crushing vessel.
- Put the dedicated crusher into the crushing vessel and close the lid tightly.
- Submerge the crushing vessel in liquid nitrogen completely with the adjustable plier etc.
- Close the lid of styrofoam container and wait until the liquid nitrogen boiled.
- Take out the crushing vessel with *3 adjustable plier etc. when the boiling settles down.
- Place the crushing vessel that taken out and take it by gloved hand.
- Place the crushing vessel on the unit and put the lid of mounting rack.
- Secure the lid of rack with the black-colored knob.
- Shaken at prescribed speed and time.
- After the shaking is complete open the lid and check inside it. *4
- Completed if the sample is crushed satisfactory. *5

*1. Desirable to use the minimum-sized polystyrene foam container that the required number crushing vessels can be immersed to minimize the amount used of liquid nitrogen.
 *2. The processing capacity of the crushing vessel is 1.2g per 1pc while better to make it to 0.5g per 1pc for plastic samples (Up to 1g polystyrene easily that can be crushed by freezing).
 *3.Wait for at least 2 minutes after the boiling is settled out to freeze the sample in the crushing vessel sufficiently.
 *4.The crushed sample might stick to the crusher so tap it with the inner wall of vessel to drop it.
 *5. If the crushing is insufficient return the crusher to the unit to freeze it again and shaken.