

# Reciprocal shaker NR-1, Rotary shaker NR-2, Double shaker NR-3

**Compact sized easy to use for bench tops and inside Incubators. 3 models: Reciprocal, Orbital, Switchable Reciprocal shaking. Options are compatible with those of Bioshaker (Small sized).**

"Bioshaker BR-21/22/23 with temp. control" series --> P.018 OD-Monitors measuring turbidity while shaking --> P.034

Other platforms can be used with each model in combination other than the combination shown below.

NR-1 with Optional shaking platform



NR-2 with Optional shaking platform



NR-3 with Optional shaking platform




## Features

- Reciprocal shaking of Test tubes and Shake flasks [NR-1]
- Orbital shaking of mainly Erlenmeyer flasks [NR-2]
- Switchable Reciprocal/Orbital shaking [NR-3]

### Compatibility of Reciprocal shaking and Vessels.

Reciprocal shaking is suitable for Test tubes and Shake flasks. In shaking culture with Shake flasks, adjust the flow rate to 1/2 of the capacity so that the liquid can be moved properly.




## Applications

- Shaking culture of Microbes inside Incubators etc.
- Shaking culture of Algae etc. inside Incubator with LED irradiation
- Staining/Decolorization of Gels and Membranes [NR-80]

### Compatibility of Orbital shaking and Vessels.


Orbital shaking is suitable for Erlenmeyer flasks. In shaking culture with Erlenmeyer flasks, adjust the flow rate to 1/5 to 1/3 of the capacity so that the liquid can be moved properly.



Model	NR-1	NR-2	NR-3
Shaking motion	Reciprocal shaking	Orbital shaking	Switchable Reciprocal/Orbital shaking
Shaking speed	20 to 200r/min		
Shaking width	10 to 40mm, Stepless (Default: 25mm)		
Platform dimensions and Capacity	290 x 250mm Erlenmeyer flask 500mL x 4, 3L x 1		
Maximum load	Approx. 5kg (Shaking platform and vessel holders included)		
Shaking duration integrator	0.1h to 999.9h (Automatic resetting function)		
Safety devices/functions	Motor protection circuit at overload, Fuse		
Dimensions / Weight	355 x 355 x 160mm, Approx. 13kg		
Power supply	AC220 to 240V/0.5A		

(\*Shaking platform and vessel holder are option.

### Dimensions (Each model common)



### Related products: OD-Monitor series and Small-sized incubator



Can be combined with the "OD-Monitor series (page 34) that can automatically measure the turbidity of culture liquid (OD 600) without stopping shaking culture. Prepare the universal shaking platform and clamps when combining.




Can be put on the lowest shelf of "In vitro Box iB-130 (page 138) for use. \* No adverse effects from vibration verified. This is an energy-saving small size incubator that can control the temp. below RT. Very useful if no incubation room when culturing or incubation that requires low temp.

Selection guide  
 Constant temperature incubator shaker  
 OD-Monitor  
 CO<sub>2</sub> incubator shaker  
 Shaker  
 Mixer  
 Rotator  
 Stirrer  
 Bead beater  
 homogenizer  
 Ultrasonic homogenizer  
 Aluminum block bath  
 Minimize bath  
 Water bath  
 Shaking water bath  
 Immersion cooler  
 Hybridation oven  
 Constant temperature  
 Centrifugal concentrator  
 Cold trap  
 Freeze dryer  
 Electrophoresis and blotting apparatus  
 Constant temperature water circulating system [Chiller]  
 Appendix

# Optional accessories and Capacity of Vessels

Small size incubator "in vitro box iB-130"-->P.138

## Spring net shaking platform and Capacity of Vessels.


Model	MR-2925			
	Vessels can be easily mounted on regardless its shapes. It can be used by removing some springs and changing the height of spring net as necessary. Universal shaking platform (w/o Clamp *It is option) recommended if vessel should be tightly held.			
Spring net Pitch	20mm			
Spring net Pitch	1 tier (Height is changeable.)			
Dimension	290x250x108Hmm			
Weight	Apporx.1.1kg			
Capacity of Vessel	φ16 mm test tubes (*1)	50 (*2)	Erlenmeyer flasks 250mL	5
	Disposable Centrifuge tubes for 50 mL	20 (*2)	Erlenmeyer flasks 300mL	5
	Erlenmeyer flasks 50mL	20	Erlenmeyer flasks 500mL	4
	Erlenmeyer flasks 100mL	10	Erlenmeyer flasks 1L	2
	Erlenmeyer flasks 200mL	9	Erlenmeyer flask 2L	1

(\*1)When below φ16mm test tubes, fix the test tube stand (incl. φ16mm) whole with the spring net.


(\*2)The number is when Angle rack tilted 45°.

(\*3)Be careful when using Shake flask as it is fragile.

## Universal shaking platform and Capacity of Flasks


Model	MT-2925			
	It is used together Clamps, Angle racks and "OD-Monitor" series (Clamps of Erlenmeyer flask and Shake flask are available). It can hold flasks tightly than Spring shaking platform.			
Dimension	290x250x23Hmm			
Weight	Apporx.0.9kg			
Capacity of Vessels	Erlenmeyer flasks 50mL	20	Erlenmeyer flasks 500mL	4
	Erlenmeyer flasks 100mL	10	Erlenmeyer flasks 1L	2
	Erlenmeyer flasks 200mL	9	Erlenmeyer flask 2L	1
	Erlenmeyer flasks 250mL	5	Erlenmeyer flask 3L	1
	Erlenmeyer flasks 300mL	5	Sakaguchi flasks 500mL	4

## Angle rack for Test tube (Mounted on MT-2925 for use.)


Model	AT-1250	AT-1620	AT-1650	AT-1820
	Capacity			
	φ12mm x50	φ16.5mm x20	φ16.5mm x50	φ18mm x20
	Qty of mounting on MT-2925			
	2 (*)	2 (*)	2 (*)	2 (*)
	Dimension (mm)			
	250x100x95H	285x70x80H	285x110x100H	320x70x80H
	Weight			
	Apporx.0.8kg	Apporx.0.6kg	Apporx.0.9kg	Apporx.0.7kg

(\*The qty when Angle rack with 160mm test tube tilted 45°.

## Angle rack for Disposable Centrifuge tubes (Mounted on MT-2925 for use.)

Model	AT-3518	
	Capacity	
	50mLx18	
	Qty of mounting on MT-2925	
	2 (When Angle rack tilted 45°)	
	Dimension/Weight	
	285x110x102Hmm, Apporx.0.8kg	
	Remarks	
	Be sure use below 50°C.	

## Stick sheet shaking platform

Model	SR-2925
	Flat bottoms such as Erlenmeyer flasks and Petri dishes can be held simply by placing them on the platform. The decreasing an adhesive force by dirt can be recovered by washing with water and drying it. Even Erlenmeyer flasks can be used up to 200r/min by using Optional Clamp belt.
Max. shaking speed	100r/min (200r/min when using Clamp belt)
Dimension/Weight	290x250x30Hmm, Apporx.1.1kg
Standard accessories	Stick sheetx1

## Stick sheet, Clamp belts



Stick sheet (ST-2925) above also available as sold separately.

Clamp belt (CB-0205/1000) is for securely holding Erlenmeyer flask on Stick sheet shaking platform.

Model	ST-2925	Model	CB-0205	CB-1000
Dimension	290x250mm	Volume of Erlenmeyer flasks	200 to 500mL	1 to 3L (*)


(\*Up to 2L of Erlenmeyer flask can be used with the NR-1/2/3 series.

## Clamps (Mounted on MT-2925 and AF-2070 for use.)

Shape	Vessels / Capacity	Model	
	Erlenmeyer flasks	50mL	<b>CF-0050</b>
		100mL	<b>CF-0100</b>
		200mL	<b>CF-0200</b>
		250mL	<b>CF-0250</b>
		300mL	<b>CF-0300</b>
		500mL	<b>CF-0500</b>
		1L	<b>CF-1000</b>
		2L	<b>CF-2000</b>
		3L	<b>CF-3000</b>
			Sakaguchi flasks


## Angle rack for Erlenmeyer flasks

(Mounted on MT-2925 \*Clamps required for use)

Model	AF-2070	AF-2570
	Capacity	
	50mL:4,100mL:3pcs, 200 to 300mL:2pcs, 500mL:1pc	50mL:5pcs, 100mL:4pcs, 200mL:3pcs, 250 to 500mL:2pcs
	Qty of mounting on MT-2925	
	2 pcs (*1)	
	Dimension	
	250x110x50Hmm	
	Weight	
	Apporx.0.8kg	

(\*1)When Angle rack tilted 30°

## Tube holders (Mounted on MT-2925 for use.)

Model	M-0150	
	Capacity	
	1.5/2.0 mL Microtubes and Micro Vials, 3 sets x 8 racks	
	Qty of mounting on MT-2925	
	3 sets x 8 racks (*1)	
	Dimension/Weight	
	285x110x102Hmm, Apporx.0.8kg	
	Remarks	
	Be sure use below 50°C.	

(\*1) Up to 3 sets can be mounted on MT-2925 in NR-1/2/3 series. The vessel tips will protrude from MT-2925 when 3pcs are attached. Thus, be sure not to make them hit the surroundings during shaking.